



## Cholera toxin BSA conjugated

---

The hapten-carrier coupled antigen paired with mouse monoclonal antibody has been validated for Enzyme Linked Immunosorbent Assay (ELISA) & Lateral Flow Assay (LF).

Catalog No.	821158
Size	
Product Category	Small Molecule Antigen
Classification	Mycotoxin
Purity	>90% (by SDS-PAGE gel)
Storage Buffer	PBS pH7.112
Applications	ELISA & Lateral Flow (LF)
Storage/Stability	-20°C / 1 year
Shipping	Gel Packs