

SABC (Rabbit IgG-AP) kit

Streptavidin-Biotin Complex (SABC) kit is specially designed for displaying the distribution of antigens on tissues and cells in immunochemistry and other immunodetection analyses. Streptavidin is a 47,000 Dalton protein purified from the bacterium Streptomyces avidinii. Streptavidin has extraordinarily strong affinity to biotin molecules. The dissociation constant (Kd) of the biotin-Streptavidin complex is on the order of ~ 10 -15 mol/L, a million times higher than the typical affinity between antigens and their antibodies. Streptavidin has very low non-specific binding to tissues and cells, due to its nearly neutral isoelectric point (pI=6.0 \sim 6.5). Therefore, immunohistochemical analyses based on Streptavidin-biotin complex has extremely low background. Furthermore, this kit has high sensitivity because each complex it generates has a large number of alkaline phosphatase and Streptavidin molecules. In brief, SABC offers high specificity, low background and ease-of-use.

Note: Rabbit IgG refers to the host of the primary antibody, not the origin of the specimen. This kit should be used on primary antibodies from Rabbit; the compnents supplied in this kits are sufficient for 100-200 assay.

| Catalog No. | 504317 |
|-------------------|--|
| Size | Kit |
| Product Category | Immuno Detection |
| Kit Components | 1) Biotinylated Goat anti-Rabbit IgG: 100uL 2) AP-Streptavidin concentrate: 100uL 3) BSA Blocking Buffer (5%): 10mL 4) Dilution buffer: 30mL 5) 1x AP Buffer: 30mL 6)BCIP [25x]: 1mL 7) NBT [25x]: 1mL |
| Storage/Stability | -20°C / 1year |