



ELISA Coating Buffer [10x]

The 10x stock solution of buffer can be diluted 10 times with deionized water to form a ready-to-use working solution of pH 9.6. While using, take an appropriate amount of coating buffer to dilute the antigen or antibody solution to the corresponding coating concentration, add 100uL to each reaction well of the ELISA plate, and coat overnight at 4°C.

Catalog No.	505017
Size	100mL / 500mL
Product Category	Immuno Reagent
Storage/Stability	Ambient/1 year
Shipping	Ambient