

Human Mothers Against Decapentaplegic Homolog 2 (Smad2/MADH2) Rapid ELISA kit

This Rapid ELISA kit applies Sandwich enzyme immunoassay technique for the quantitative detection. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human Smad2/MADH2. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human Smad2/MADH2 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human Smad2/MADH2, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Human Smad2/MADH2. You can calculate the concentration of Human Smad2/MADH2 in the samples by comparing the OD of the samples to the standard curve.

Catalog No.	3522706	
Size	96-Well Tests	
Product Category	Rapid ELISA	
Reactivity	Human	
Sample	Serum; Plasma; Tissue Homogenate; Other Biological Samples	
Assay Method	Sandwich ELISA	
Assay Duration	1.5H	
Sensitivity	46.88 pg/mL	
Standard Curve Range	78.13-5000 pg/mL	

C+	/C+-	La Hilland
Stora	ge/Sta	IDIIITY

-20°C/1 year; 4°C/6 months

www.realgenelabs.com

For Research Use only