



Rat Cutaneous T-Cell Attracting Chemokine (CTACK) Rapid ELISA kit

This Rapid ELISA kit applies Sandwich enzyme immunoassay technique for the quantitative detection. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat CTACK. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Rat CTACK and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Rat CTACK, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of $450 \text{ nm} \pm 2 \text{ nm}$. The OD value is proportional to the concentration of Rat CTACK. You can calculate the concentration of Rat CTACK in the samples by comparing the OD of the samples to the standard curve.

Catalog No.	3541017
Size	96-Well Tests
Product Category	Rapid ELISA
Reactivity	Rat
Sample	Serum; Plasma; Tissue Homogenate; Other Biological Samples
Assay Method	Sandwich ELISA
Assay Duration	1.5H
Sensitivity	37.5 pg/mL
Standard Curve Range	62.5-4000 pg/mL
Storage/Stability	-20°C/1 year; 4°C/6 months

