



## REAL-Detect Plant PAT gene qPCR kit

---

REAL-Detect Plant PAT gene qPCR kit applies qualitative fluorescence PCR method for gene detection. The kit is supplied with reaction mix, positive control, negative control and re-dissolved diluent. The components supplied in the kit are in lyophilized format (stable at 25-30°C for 2 weeks).

Kit is compatible with all the standard Real-Time PCR systems, including:

ABI: 7000, 7300, 7500, 7900, Step One, Step One Plus, QuantStudio

Agilent: Mx3000P, 3005P

Bio-Rad: CFX 96, iCycler iQ4, iCycleriQ5

Cepheid: Smart Cyclor-II

MJ: Option2, Chromo4

Qiagen: Rotor Gene 6000

Roche: LightCycler680

|                   |                      |
|-------------------|----------------------|
| Catalog No.       | 950019               |
| Size              | 48 Tests             |
| Product Category  | Transgenic Plant PCR |
| Storage/Stability | -20°C / 1 year       |
| Shipping          | Ambient              |