

## REAL-Detect Porcine Circovirus Type 2 (PCV-2) qPCR kit

REAL-Detect Porcine Circovirus Type 2 (PCV-2) qPCR kit applies qualitative fluorescence PCR method for pathogen detection. The kit is supplied with reaction mix, positive control, negative control and re-dissolved diluent. The components supplied in the kit are in lyophilized format (stable at 25-30°C for 2 weeks).

Kit is compatible with all the standard Real-Time PCR systems, including: ABI: 7000, 7300, 7500, 7900, Step One, Step One Plus, QuantStudio Agilent: Mx3000P, 3005P Bio-Rad: CFX 96, iCycler iQ4, iCycleriQ5 Cepheid: Smart Cycler-II MJ: Option2, Chromo4 Qiagen: Rotor Gene 6000

Roche: LightCycler680

Catalog No.	921117
Size	48 Tests
Product Category	Veterinary PCR
Storage/Stability	-20°C / 1 year
Shipping	Ambient
Pack Size / Price	Inquire

www.realgenelabs.com

For Research Use only