

Rapid Stain-free Gel Preparation kit [6%]

Rapid Stain-free Gel Preparation kit is suitable for the rapid preparation of polyacrylamide gels. This product adopts the pre-mixed formula; only needs to mix the reagents two by two, and then add the coagulant to gel. After pouring the lower layer of glue, there is no need for a liquid seal, and the upper layer of glue can be poured directly, which is convenient and quick. The matching top layer tape is colored, which is convenient for spotting and distinguishing different gels. The gel contains fluorescent dyes, which can directly observe the protein separation effect and membrane transfer efficiency in the imager. The improved coagulant provided has good stability and coagulation acceleration effect, and there is no need to add TEMED during the compounding process. The gel made by This kit can be used for denaturing PAGE gel electrophoresis.

Catalog No.	502316
Size	20 preps / 50 preps
Product Category	SDS-PAGE Buffer & Reagent
Storage/Stability	2 ~ 8°C/1 year
Shipping	Gel Packs

www.realgenelabs.com

For Research Use only