

## Rapid Blood Protein Extraction Kit (spin-column method)

Blood Protein rapid extraction kit (spin-column method) is suitable for extracting total protein from all kinds of anticoagulant or non-anticoagulant whole blood samples. The extraction process is simple and convenient, and can be completed within 1 hour. The kit contains a protease inhibitor mixture that prevents protease from degrading the protein and provides a guarantee for extracting high purity proteins. This kit contains a unique formulation that effectively dissolves cell membrane components, including the plasma membrane, the nuclear membrane, and various organelle membranes. The kit contains a protease inhibitor mixture that prevents protease from degrading the protein and provides a guarantee for extracting high purity proteins. The protein extracted from this kit can be used for western-blotting, protein electrophoresis, immunoprecipitation, ELISA, transcriptional activity analysis, Gelshift gel retardation test, enzyme activity determination and other downstream protein research experiments. The protein extracted by this kit is an active protein with natural protein conformation. This kit does not contain EDTA and is compatible with downstream applications such as metal chelation and chromatography.

Catalog No.	501052
Size	50 preps
Product Category	Protein Extraction
Storage/Stability	2-8°C/lyear
Shipping	Ambient

www.realgenelabs.com

For Research Use only